

V2 Collection Vessel Close Down and Completion

Chart 4.

Completing the test and sealing the sample

1. Removing scope from collection vessel

With insertion tube still inserted into the silicon insert, carefully disengage the silicon insert with scope attached.



2. Removal and engagement of lid to seal sample

Gently remove the collection set from the aluminium upstand
Place on a flat surface

With the palm of your hand and without tilting the sample push down on the top to engage the jar with the lid

If the jar is off centre to the lid, gently wiggle till you feel it engage, twist and lock the lid onto the jar.



3. Removal of sample from vessel

Once the lid has engaged and sample sealed, place your hands around the silicon base and jar and with a twisting wiggle motion start to pull the sample from the vessel. You can lubricate the jar a little with your alcohol towel if it helps.

Vessel and silicon fittings can be discarded or recycled.



4. Removing scope from silicon insert

At end of test - Hold silicon insert, remove scope tip, discard silicon insert.
Utilise alcohol wipes as lubricant if required around IT, wiping down from end of silicon insert to tip of scope as you remove.

Label and bag up sample ready for the lab.

