

ENDOTEST V2 MICROTESTING FLOW CHART

TEST SAMPLE ORDER FOR TESTING DUAL CHANNEL SCOPE USING HAND PUMP

Chart 3.

A. Wetting – Filling biopsy and suction channels

- All lines should be attached to the Endoscope as per accessory kit instructions Chart 2.
- Close all clamps

Open Yellow Clamp

Using hand pump, ½ pump handle and as soon as fluid appears in the vessel. CLOSE CLAMP. Your channel is now filled with fluid ready for brushing and collection.

Be careful not to over pump and close clamps as soon as you witness fluid enter the collection vessel
Too much pressure and not closing the clamps will fill your collection vessel too quick.

Open 1st White Clamp

Using hand pump ½ pump, again as soon as you witness fluid enter vessel.
Close White clamp

Open 2nd White Clamp

Using hand pump ½ pump, again as soon as you witness fluid enter vessel.
Close White clamp

B. Collection of sample from air and water Channels

Open Red Clamp

Using hand pump, pump handle until you witness a good stream of fluid enter collection vessel 10ml. Can measure using side of vessel.
Close Red clamp.

Open Blue Clamp

Using hand pump, pump handle until you witness a good stream of fluid enter collection vessel 10ml. Can measure using side of vessel.
Close Blue clamp.

C. Brushing (prior to entering an aseptic method)

First – check your scope tip, angulate down and lock to ensure brush can enter the pot with ease.

- Remove suction cylinder cap and biopsy port connectors, place on aseptic field
- Prepare your brush and gloves for easy access.
- Wash hands and put on sterile gloves.
(following the aseptic method)
- Obtain brush, touch brush and nothing else then proceed to brush all brush able channels
 - Suction channel, push through until brush is sighted in the silicon tube connected at proximal connector - Remove
 - Suction link channel, enter suction port at 45° push through to approx. 200mm past biopsy port - Remove
 - Biopsy channel, enter biopsy port and push through until the brush is sighted in the fluid in the collection pot – wash in fluid.
 - Repeat for 2nd biopsy channel
 - Remove and discard brush.
 - Replace suction cylinder cap and biopsy port connector.

Re angulate your scope tip so it can be removed easily at the end of the test

D. Final collection of brushed sample

Open Yellow clamp

Using hand pump, pump handle to pull brushed sample and empty fluid from channel.
Remove the bottle line from water bottle to avoid over filling
You should hear the water being pulled from the channel
Once you witness fluid being pulled from the channel. Close Yellow Clamp

Open White clamps

Pump if required – this allows remaining fluid to be drained from the channels
You should hear the water being pulled from the channels
Test is finished.
Follow instructions Chart 4. to seal the vessel