



## ENDOTEST V2 MICROBIOLOGICAL SCOPE COMPITABILITY LIST

**July 2025 - Jun 2026**

*If you cannot find your scope model listed please contact us for advice*



**\*\* Some 3.5-4.1 IT sized that used a 1.4mm kit codes have changed to 4-8mm due to new silicon inserts  
Please check codes before ordering in case of the change**

**OLYMPUS**

|                                                                                                                                     | Ext Diam      | Channel Diam/WL-TL         | Collection Pot Size | Brush Size | Endoworx V2 Order Code |
|-------------------------------------------------------------------------------------------------------------------------------------|---------------|----------------------------|---------------------|------------|------------------------|
| <b>#L Denotes supplied with a 5mm diameter brush for use on channels &gt; 2.6mm. #R 3.5-4.1mm IT with 1.5-1.7 Channel 4-8mm kit</b> |               |                            |                     |            |                        |
| <b>Bronchoscopes</b>                                                                                                                |               |                            |                     |            |                        |
| BF-H1100                                                                                                                            | 4.9mm         | 2.2mm/600                  | (4-8mm)             | 3mm 1500   | 727-ET2405             |
| BF-H190                                                                                                                             | 5.5mm         | 2.0mm/600-890              | (4-8mm)             | 3mm 1500   | 727-ET2405             |
| BF-IT160/IT180                                                                                                                      | 6.0mm         | 3.0mm/600-870              | (4-8mm)             | 5mm 2400   | 727-ET2405L            |
| BF-IT180Q                                                                                                                           | 6.0mm         | 3.0mm/600-870              | (4-8mm)             | 5mm 2400   | 727-ET2405L            |
| BF-ITH1100                                                                                                                          | 5.8mm         | 3.0mm/600                  | (4-8mm)             | 5mm 2400   | 727-ET2405L            |
| BF-ITH190                                                                                                                           | 6.2mm         | 2.8mm/600-860              | (4-8mm)             | 5mm 2400   | 727-ET2405L            |
| BF-MP190F                                                                                                                           | 3.7mm         | 1.7mm/600                  | (4-8mm)             | 1.5mm 1500 | 727-ET2405R            |
| BF-P180                                                                                                                             | 4.8mm         | 2.0mm/600-870              | (4-8mm)             | 3mm 1500   | 727-ET2405             |
| July 2025 - Jun 2026                                                                                                                | 4.2mm         | 2.0mm/600-870              | (4-8mm)             | 3mm 1500   | 727-ET2405             |
| BF-PE2                                                                                                                              | 5.0mm         | 2.2mm/600-890              | (4-8mm)             | 3mm 1500   | 727-ET2405             |
| BF-Q180                                                                                                                             | 5.5mm         | 2.0mm/600-870              | (4-8mm)             | 3mm 1500   | 727-ET2405             |
| BF-Q190                                                                                                                             | 4.9mm         | 2.0mm/600-870              | (4-8mm)             | 3mm 1500   | 727-ET2405             |
| BF-TE2                                                                                                                              | 6.0mm         | 2.8mm/600-890              | (4-8mm)             | 5mm 2400   | 727-ET2405L            |
| BF-UC180F/190F                                                                                                                      | 6.9mm         | 2.2mm/600-890              | (4-8mm)             | 3mm 1500   | 727-ET2403             |
| BF-XP60                                                                                                                             | 2.8mm         | 1.2mm/600-950              | (1-4mm)             | 1.5mm 1500 | 727-ET2105             |
| BF-XP160                                                                                                                            |               |                            |                     |            | 727-ET2105             |
| BF-XP190                                                                                                                            | 3.1mm         | 1.2mm/600-950              | (1-4mm)             | 1.5mm 1500 | 727-ET2105             |
| BF-XT190                                                                                                                            | 6.1mm         | 3.2mm/600-870              | (4-8mm)             | 5mm 2400   | 727-ET2405L            |
| <b>Colonoscopes</b>                                                                                                                 |               |                            |                     |            |                        |
| CF-EZ1500DL/I                                                                                                                       | 13.2mm        | 3.7mm/1330-1690, 1680-1980 | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| CF-H180AL/I                                                                                                                         | 13.9mm        | 3.7mm/1330-1690, 1680-1980 | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| CF-H180DL/I                                                                                                                         | 13.9mm        | 3.7mm/1330-1690, 1680-1980 | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| CF-H190DI/L                                                                                                                         | 13.2mm        | 3.7mm/1330-1690, 1680-1980 | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| CF-HQ190L/I                                                                                                                         | 13.2mm        | 3.7mm/1330-1690, 1680-1980 | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| CF-Q180AL/I                                                                                                                         | 13.2mm        | 3.7mm/1680-1980            | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| CF-XZ1200L/I                                                                                                                        | 13.2mm        | 3.7mm/1330-1690, 1680-1980 | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| <b>Cystoscopes</b>                                                                                                                  |               |                            |                     |            |                        |
| CYF-2                                                                                                                               | 5.3mm - 5.5mm | 2.2mm/380-660              | (4-8mm)             | 3mm 1500   | 727-ET2410             |
| CYF-4/5/VH/VAH                                                                                                                      | 5.5mm         | 2.4mm/380-700              | (4-8mm)             | 3mm 1500   | 727-ET2410             |
| CYF-4A/5A                                                                                                                           | 5.2mm         | 2.4mm/380-700              | (4-8mm)             | 3mm 1500   | 727-ET2409             |
| CYF-V2/VA2                                                                                                                          | 5.5mm         | 2.2mm/380-660              | (4-8mm)             | 3mm 1500   | 727-ET2409             |
| <b>Choledochoscopes</b>                                                                                                             |               |                            |                     |            |                        |
| CHF-CB30L                                                                                                                           | S.D 2.8mm     | 1.2mm/700-1010             | (1-4mm)             | 1.5mm 1500 | 727-ET2110             |
| CHF-CB30S                                                                                                                           | S.D 2.8mm     | 1.2mm/450-760              | (1-4mm)             | 1.5mm 1500 | 727-ET2110             |
| CHF-V                                                                                                                               | 5.2mm         | 2.0mm/380-660              | (4-8mm)             | 3mm 1500   | 727-ET2409             |
| <b>Duodenoscopes</b>                                                                                                                |               |                            |                     |            |                        |
| TJF-Q180V                                                                                                                           | 13.7mm        | 4.2mm/1240-1550            | (8-14mm)            | 5mm 2400   | 727-ET2815             |
| TJF-Q190V                                                                                                                           | 13.5mm        | 4.2mm/1240-1550            | (8-14mm)            | 5mm 2400   | 727-ET2815             |
| <b>Entroscope</b>                                                                                                                   |               |                            |                     |            |                        |
| SIF-H190                                                                                                                            | 9.2mm         | 3.2mm/2000-2280            | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| SIF-Q180                                                                                                                            | 9.2mm         | 2.8mm/2000-2345            | (8-14mm)            | 5mm 2400   | 727-ET2801             |

**OLYMPUS**

|                         |     | Ext Diam | Channel Diam/WL-TL         | Collection Pot Size | Brush Size | Endoworx V2 Order Code |
|-------------------------|-----|----------|----------------------------|---------------------|------------|------------------------|
| <b>Gastrosopes</b>      |     |          |                            |                     |            |                        |
| GIF-1TH190              |     | 10.9mm   | 3.7mm/1030-1345            | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| GIF-2TH180              |     | 12.6mm   | 2.8mm/3.7mm/1030-1350      | (8-14mm)            | 5mm 2400   | 727-ET2804             |
| GIF-EZ1500              |     | 9.6mm    | 2.8mm/1030                 | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| GIF-H180                |     | 9.8mm    | 2.8mm/1030-1345            | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| GIF-H180AL              |     | 9.8mm    | 3.7mm/1330-1680            | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| GIF-H180J               |     | 9.9mm    | 2.8mm/1030-1345            | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| GIF-H190                |     | 9.2mm    | 2.8mm/1030-1350            | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| GIF-HQ190               |     | 9.9mm    | 2.8mm/1030-1350            | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| GIF-N180                |     | 4.9mm    | 2.0mm/1100-1420            | (4-8mm)             | 5mm 2400   | 727-ET2401             |
| GIF-Q180                |     | 8.8mm    | 2.8mm/1030-1340            | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| GIF-XP160               |     | 5.9mm    | 2.0mm/1030-1345            | (4-8mm)             | 5mm 2400   | 727-ET2401             |
| GIF-XP180N              |     | 5.5mm    | 2.0mm/1100-1420            | (4-8mm)             | 5mm 2400   | 727-ET2401             |
| GIF-XP190N              |     | 5.4mm    | 2.2mm/1100-1420            | (4-8mm)             | 5mm 2400   | 727-ET2401             |
| GIF-XZ1200              |     | 9.9mm    | 2.8mm/1030-1350            | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| <b>Hysteroscope</b>     |     |          |                            |                     |            |                        |
| HYF-XP                  |     | 3.1mm    | 1.2mm/240-540              | (1-4mm)             | 1.5mm 1500 | 727-ET2110             |
| <b>Intubation</b>       |     |          |                            |                     |            |                        |
| LF-DP                   | S.D | 3.1mm    | 1.2mm/600-855              | (1-4mm)             | 1.5mm 1500 | 727-ET2111             |
| LF-GP                   | S.D | 4.1mm    | 1.5mm/600-855              | (4-8mm)             | 1.5mm 1500 | 727-ET2411R            |
| LF-TP                   |     | 5.2mm    | 2.6mm/600-855              | (4-8mm)             | 3mm 1500   | 727-ET2411             |
| LF-V                    | S.D | 4.1mm    | 1.2mm/600-855              | (4-8mm)             | 1.5mm 1500 | 727-ET2409R            |
| <b>Paediatric</b>       |     |          |                            |                     |            |                        |
| PCF-H180AI              |     | 11.3mm   | 3.2mm/1330-1690            | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| PCF-H180AL              |     | 11.7mm   | 3.2mm/1680-2000            | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| PCF-H190DL/I            |     | 11.7mm   | 3.2mm/1330-1690, 1680-2000 | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| PCF-H190L/I             |     | 11.7mm   | 3.2mm/1330-1690, 1680-2000 | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| PCF-H190TL/I            |     | 10.5mm   | 3.2mm/1330-1655, 1680-2005 | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| PCF-HQ190L/I            |     | 11.8mm   | 3.2mm/1330-1655, 1680-2005 | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| PCF-PH190L/I            |     | 9.7mm    | 3.2mm/1330-1690, 1680-2000 | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| PCF-Q160AI              |     | 11.5mm   | 3.7mm/1330-1690            | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| PCF-Q160AL              |     | 11.5mm   | 3.2mm/1680-2000            | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| PCF-Q180AI              |     | 11.3mm   | 3.2mm/1680-2000            | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| PCF-Q180AL              |     | 11.3mm   | 3.2mm/1680-2000            | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| <b>Pleuroscope</b>      |     |          |                            |                     |            |                        |
| LTF-160                 |     | 7.0mm    | 2.8mm/270-520              | (4-8mm)             | 3mm 1500   | 727-ET2405             |
| <b>Rhinolaryngo</b>     |     |          |                            |                     |            |                        |
| ENF-T3                  |     | 4-8mm    | 2.2mm/365-635              | (4-8mm)             | 3mm 1500   | 727-ET2411             |
| ENF-VT2/VT3             |     | 4.8mm    | 2.0mm/365-635              | (4-8mm)             | 3mm 1500   | 727-ET2414             |
| <b>Ultrasound/Sonic</b> |     |          |                            |                     |            |                        |
| GF-UCT180               |     | 14.6mm   | 3.7mm/1250-1575            | (8-14mm)            | 5mm 2400   | 727-ET2808             |
| GF-UE160-AL5            |     | 13.8mm   | 2.2mm/1250-1555            | (8-14mm)            | 5mm 2400   | 727-ET2807             |
| GF-UE190                |     | 13.4mm   | 2.2mm/1250-1563            | (8-14mm)            | 5mm 2400   | 727-ET2807             |
| TGF-UC180J              |     | 14.6mm   | 3.7mm/1245-1560            | (8-14mm)            | 5mm 2400   | 727-ET2801             |
| <b>Urethra scopes</b>   |     |          |                            |                     |            |                        |
| URF-P                   | S.D | 3.6mm    | 1.2mm/700-1050             | (4-8mm)             | 1.5mm 1500 | 727-ET2410R            |
| URF-P2                  | S.D | 3.3mm    | 1.2mm/700-1050             | (1-4mm)             | 1.5mm 1500 | 727-ET2110             |
| URF-P5, P6, P7          | S.D | 2.8mm    | 1.2mm/700-1050             | (1-4mm)             | 1.5mm 1500 | 727-ET2110             |
| URF-V/V3                | S.D | 3.3mm    | 1.2mm/670-980              | (1-4mm)             | 1.5mm 1500 | 727-ET2110             |

**PENTAX**

|                                                                                                                           | Ext Diam | Channel Diam/WL-TL | Collection Pot Size | Brush Size | Endoworx V2 Order Code |             |
|---------------------------------------------------------------------------------------------------------------------------|----------|--------------------|---------------------|------------|------------------------|-------------|
| #L Denotes supplied with a 5mm diameter brush for use on channels > 2.6mm. #R 3.5-4.1mm IT with 1.5-1.7 Channel 4-8mm kit |          |                    |                     |            |                        |             |
| Bronchoscopes Fibre                                                                                                       |          |                    |                     |            |                        |             |
| FB10V                                                                                                                     | S.D      | 3.5mm              | 1.15mm/600-900      | (4-8mm)    | 1.5mm 1500             | 727-ET2452R |
| FB15BS                                                                                                                    |          | 4.8mm              | 2.0mm/600-880       | (4-8mm)    | 3mm 1500               | 727-ET2452  |
| FB15RBS                                                                                                                   |          | 4.9mm              | 1.95mm/600-880      | (4-8mm)    | 3mm 1500               | 727-ET2452  |
| FB15V                                                                                                                     |          | 4.9mm              | 2.1mm/600-900       | (4-8mm)    | 3mm 1500               | 727-ET2452  |
| FB18RBS                                                                                                                   |          | 6.0mm              | 2.55mm/600-880      | (4-8mm)    | 3mm 1500               | 727-ET2452  |
| FB18V/19V                                                                                                                 |          | 6.0mm              | 2.7mm/600-900       | (4-8mm)    | 5mm 2400               | 727-ET2452L |
| FB19TV                                                                                                                    |          | 6.2mm              | 3.0mm/600-900       | (4-8mm)    | 5mm 2400               | 727-ET2452L |
| FB8V                                                                                                                      | S.D      | 2.8mm              | 1.15mm/600-900      | (1-4mm)    | 1.5mm 1500             | 727-ET2152  |
| Bronchoscopes Video                                                                                                       |          |                    |                     |            |                        |             |
| EB1170K                                                                                                                   | S.D      | 3.7mm              | 1.2mm/600-860       | (4-8mm)    | 1.5mm 1500             | 727-ET2452R |
| EB11J10                                                                                                                   |          | 3.7mm              | 1.2mm/600-860       | (4-8mm)    | 1.5mm 1500             | 727-ET2453R |
| EB15J10                                                                                                                   |          | 5.5mm              | 2.0mm/600-900       | (4-8mm)    | 3mm 1500               | 727-ET2453  |
| EB1575K                                                                                                                   |          | 5.2mm              | 2.0mm/600-860       | (4-8mm)    | 3mm 1500               | 727-ET2452  |
| EB1970K                                                                                                                   |          | 6.2mm              | 2.8mm/600-900       | (4-8mm)    | 5mm 2400               | 727-ET2452L |
| EB1970TK                                                                                                                  |          | 6.3mm              | 3.2mm/600-860       | (4-8mm)    | 5mm 2400               | 727-ET2452L |
| EB1975K                                                                                                                   |          | 6.4mm              | 2.8mm/600-860       | (4-8mm)    | 5mm 2400               | 727-ET2452L |
| EB1990i                                                                                                                   |          | 6.4mm              | 1.2mm/600-860       | (4-8mm)    | 1.5mm 1500             | 727-ET2452R |
| EB19J10                                                                                                                   |          | 6.4mm              | 2.8mm/600-900       | (4-8mm)    | 5mm 2400               | 727-ET2453L |
| Bronchoscopes Video Ultrasound                                                                                            |          |                    |                     |            |                        |             |
| EB1970UK                                                                                                                  |          | 6.3mm              | 2.0mm/600-860       | (4-8mm)    | 3mm 1500               | 727-ET2455  |
| Colonoscopes                                                                                                              |          |                    |                     |            |                        |             |
| EC2990Fi/Li/Mi                                                                                                            |          | 10.2mm             | 2.8mm/1700-2000     | (8-14mm)   | 5mm 2400               | 727-ET2850  |
| EC3490Fi/Li/Mi/TFi/TLi/TMi                                                                                                |          | 11.6mm             | 3.2mm/1300-2023     | (8-14mm)   | 5mm 2400               | 727-ET2850  |
| EC3490FK/LK                                                                                                               |          | 11.6mm             | 3.8mm/1500-2023     | (8-14mm)   | 5mm 2400               | 727-ET2850  |
| EC34-i10F/i10L/i10M                                                                                                       |          | 11.6mm             | 3.8mm/1300-2000     | (8-14mm)   | 5mm 2400               | 727-ET2850  |
| EC34-i10TF/i10TL/i10TM                                                                                                    |          | 11.6mm             | 3.2mm/1300-2000     | (8-14mm)   | 5mm 2400               | 727-ET2850  |
| EC3890Fi/Li/Mi                                                                                                            |          | 13.2mm             | 3.8mm/1300-2023     | (8-14mm)   | 5mm 2400               | 727-ET2850  |
| EC3890FK/FK2/K/MK                                                                                                         |          | 13.2mm             | 4.2mm/1300-2023     | (8-14mm)   | 5mm 2400               | 727-ET2850  |
| EC3890FZi/LZi/MZi                                                                                                         |          | 13.2mm             | 3.8mm/1300-2025     | (8-14mm)   | 5mm 2400               | 727-ET2850  |
| EC38-i10F/i10L/i10M                                                                                                       |          | 13.2mm             | 3.8mm/1300-2000     | (8-14mm)   | 5mm 2400               | 727-ET2850  |
| EC38-i10F2/i10M2                                                                                                          |          | 13.2mm             | 3.8mm/1300-2000     | (8-14mm)   | 5mm 2400               | 727-ET2850  |
| Fiberoptic Cholangioscope                                                                                                 |          |                    |                     |            |                        |             |
| FCP-9P                                                                                                                    |          | 3.1mm              | 1.15mm/1900-2100    | (1-4mm)    | Brush not supplied     | 727-ET2152  |
| Cystoscopes                                                                                                               |          |                    |                     |            |                        |             |
| FCY-15RBS                                                                                                                 |          | 4.9mm              | 2.1mm/400-680       | (4-8mm)    | 3mm 1500               | 727-ET2459  |
| ECY-1570/1575K                                                                                                            |          | 5.2mm              | 2.0mm/400-660       | (4-8mm)    | 3mm 1500               | 727-ET2459  |
| Duodenoscopes                                                                                                             |          |                    |                     |            |                        |             |
| ED34-i10T/2                                                                                                               |          | 13.0mm             | 4.2mm/1250-1550     | (8-14mm)   | 5mm 2400               | 727-ET2856  |
| ED3470TK / ED3490TK                                                                                                       |          | 13.2mm             | 4.2mm/1250-1550     | (8-14mm)   | 5mm 2400               | 727-ET2856  |
| Gastrosopes                                                                                                               |          |                    |                     |            |                        |             |
| EG1690K                                                                                                                   |          | 5.4mm              | 2.0mm/1100-1423     | (4-8mm)    | 3mm 1500               | 727-ET2454  |
| EG16-K10/EG17-J10                                                                                                         |          | 5.4mm              | 2.0mm/1100-1416     | (4-8mm)    | 3mm 1500               | 727-ET2454  |
| EG1870K                                                                                                                   |          | 6.0mm              | 2.0mm/1050-1400     | (4-8mm)    | 3mm 1500               | 727-ET2450  |
| EG2490K                                                                                                                   |          | 8.0mm              | 2.4mm/1050-1373     | (4-8mm)    | 3mm 1500               | 727-ET2454  |
| EG2790i/2990i                                                                                                             |          | 9.0mm              | 2.8mm/1050-1373     | (8-14mm)   | 5mm 2400               | 727-ET2850  |
| EG2790K                                                                                                                   |          | 9.0mm              | 2.8mm/1050-1373     | (8-14mm)   | 5mm 2400               | 727-ET2850  |
| EG27-i10                                                                                                                  |          | 9.0mm              | 2.8mm/1050-1366     | (8-14mm)   | 5mm 2400               | 727-ET2850  |
| EG2990K                                                                                                                   |          | 9.8mm              | 2.8mm/1050-1373     | (8-14mm)   | 5mm 2400               | 727-ET2850  |
| EG2990Zi                                                                                                                  |          | 9.8mm              | 2.8mm/1050-1400     | (8-14mm)   | 5mm 2400               | 727-ET2850  |
| EG29-i10/i10N                                                                                                             |          | 9.8mm              | 3.2mm/1050-1366     | (8-14mm)   | 5mm 2400               | 727-ET2850  |
| EG3490K                                                                                                                   |          | 11.6mm             | 3.8mm/1050-1373     | (8-14mm)   | 5mm 2400               | 727-ET2850  |
| EG34-i10                                                                                                                  |          | 11.00mm            | 3.8mm/1050-1400     | (8-14mm)   | 5mm 2400               | 727-ET2850  |
| FG29V                                                                                                                     |          | 9.8mm              | 2.8mm/1050-1395     | (8-14mm)   | 5mm 2400               | 727-ET2850  |

**PENTAX**

|                          |     | Ext Diam | Channel Diam/WL-TL  | Collection Pot Size | Brush Size | Endoworx V2 Order Code |
|--------------------------|-----|----------|---------------------|---------------------|------------|------------------------|
| <b>Dual Channel</b>      |     |          |                     |                     |            |                        |
| EG3890TK                 |     | 12.8mm   | 3.8,2.8mm/1050-1383 | (8-14mm)            | 5mm 2400   | 727-ET2851             |
| EC3890TLK/TFK            |     | 13.2mm   | 3.8,2.8mm/1700-2033 | (8-14mm)            | 5mm 2400   | 727-ET2851             |
| <b>Intubation Scopes</b> |     |          |                     |                     |            |                        |
| Fi-10RBS                 | S.D | 3.5mm    | 1.3mm/600-900       | (4-8mm)             | 1.5mm 1500 | 727-ET2452R            |
| Fi-13RBS                 | S.D | 4.1mm    | 1.7mm/600-900       | (4-8mm)             | 1.5mm 1500 | 727-ET2452R            |
| Fi-16RBS                 |     | 5.2mm    | 2.5mm/600-900       | (4-8mm)             | 3mm 1500   | 727-ET2452             |
| <b>Rhinolaryngoscope</b> |     |          |                     |                     |            |                        |
| FNL-15RP3                |     | 4.9mm    | 2.1mm/300-580       | (4-8mm)             | 3mm 1500   | 727-ET2452             |
| VNL-1570STK              |     | 4.9mm    | 2.00mm/300-560      | (4-8mm)             | 3mm 1500   | 727-ET2452             |
| <b>Ultrasound (EUS)</b>  |     |          |                     |                     |            |                        |
| EG3270UK Linear          |     | 10.8mm   | 2.8mm/1250-1560     | (8-14mm)            | 5mm 2400   | 727-ET2858             |
| EG34-J10U Linear         |     | 14.3mm   | 2.8mm/1250-1566     | (8-14mm)            | 5mm 2400   | 727-ET2858             |
| EG3670URK Rad            |     | 12.1mm   | 2.4mm/1250-1560     | (8-14mm)            | 5mm 2400   | 727-ET2857             |
| EG36-J10UR Rad           |     | 10.4mm   | 2.4mm/1250-1566     | (8-14mm)            | 5mm 2400   | 727-ET2857             |
| EG3870UTK Linear         |     | 12.8mm   | 3.8mm/1250-1560     | (8-14mm)            | 5mm 2400   | 727-ET2858             |
| EG38-J10UT Linear        |     | 12.9mm   | 4.00mm/1250-1566    | (8-14mm)            | 5mm 2400   | 727-ET2858             |

OTHER

| Fujinon                                                                                                                   | Ext Diam | Channel Diam/WL-TL | Collection Pot Size | Brush Size | Endoworx V2 Order Code |            |
|---------------------------------------------------------------------------------------------------------------------------|----------|--------------------|---------------------|------------|------------------------|------------|
| #L Denotes supplied with a 5mm diameter brush for use on channels > 2.6mm. #R 3.5-4.1mm IT with 1.5-1.7 Channel 4-8mm kit |          |                    |                     |            |                        |            |
| Balloon Scopes                                                                                                            |          |                    |                     |            |                        |            |
| EI-580BT                                                                                                                  | 9.4mm    | 3.2mm/1550-1850    | (8-14mm)            | 5mm 2400   | 727-ET2887             |            |
| EN-450T5/W                                                                                                                | 9.4mm    | 2.8mm/2000-2300    | (8-14mm)            | 5mm 2400   | 727-ET2883             |            |
| EN-580T                                                                                                                   | 9.4mm    | 3.2mm/2000-2300    | (8-14mm)            | 5mm 2400   | 727-ET2887             |            |
| EN-580XP                                                                                                                  | 7.5mm    | 2.2mm/2000-2300    | (4-8mm)             | 5mm 2400   | 727-ET2487             |            |
| Bronchoscopes                                                                                                             |          |                    |                     |            |                        |            |
| EB-530H                                                                                                                   | 5.4mm    | 2.0mm/600-870      | (4-8mm)             | 3mm 1500   | 727-ET2493             |            |
| EB-530P                                                                                                                   | 3.8mm    | 1.2mm/600-890      | (4-8mm)             | 1.5mm 1500 | 727-ET2493R            |            |
| EB-530S                                                                                                                   | 4.9mm    | 2.0mm/600-870      | (4-8mm)             | 3mm 1500   | 727-ET2493             |            |
| EB-530T                                                                                                                   | 5.9mm    | 2.8mm/600-870      | (4-8mm)             | 3mm 1500   | 727-ET2493             |            |
| July 2025 - Jun 2026                                                                                                      | 6.2mm    | 3.2mm/600-870      | (4-8mm)             | 3mm 1500   | 727-ET2493             |            |
| EB-580S                                                                                                                   | 5.3mm    | 2.2mm/600-870      | (4-8mm)             | 3mm 1500   | 727-ET2493             |            |
| EB-580T                                                                                                                   | 5.9mm    | 2.8mm/600-870      | (4-8mm)             | 3mm 1500   | 727-ET2493             |            |
| EB-710P                                                                                                                   | 4.1mm    | 2.0mm/600-880      | (4-8mm)             | 3mm 1500   | 727-ET2493             |            |
| Colonoscopes                                                                                                              |          |                    |                     |            |                        |            |
| Scopes that have the jet channel on the control head require the use of a jet channel adapter                             |          |                    |                     |            |                        |            |
| EC-530FM/FI/FL                                                                                                            | 12.80mm  | 3.8mm/1520-1990    | (8-14mm)            | 5mm 2400   | 727-ET2890             |            |
| EC-530IT/LT/MT                                                                                                            | 12.8mm   | 4.2mm/1330-1990    | (8-14mm)            | 5mm 2400   | 727-ET2890             |            |
| EC-530MP/LP                                                                                                               | 11.00mm  | 3.2mm/1330 to 1990 | (8-14mm)            | 5mm 2400   | 727-ET2890             |            |
| EC-530WI3/WL3/WM3                                                                                                         | 12.80mm  | 3.8mm/1330 to 1990 | (8-14mm)            | 5mm 2400   | 727-ET2890             |            |
| EC-580RDM/L                                                                                                               | 10.50mm  | 3.2mm/1330-1990    | (8-14mm)            | 5mm 2400   | 727-ET2890             |            |
| EC-590WI4/WL4/WM4                                                                                                         | 12.80mm  | 3.8mm/1330-1990    | (8-14mm)            | 5mm 2400   | 727-ET2890             |            |
| EC-590ZW3L/ZW3M                                                                                                           | 12.80mm  | 3.8mm/1330-1990    | (8-14mm)            | 5mm 2400   | 727-ET2890             |            |
| EC-600WI/WL/WM                                                                                                            | 12.00mm  | 3.8mm/1330 to 1990 | (8-14mm)            | 5mm 2400   | 727-ET2890             |            |
| EC-600ZWL/M                                                                                                               | 12.80mm  | 3.8mm/1330 to 1990 | (8-14mm)            | 5mm 2400   | 727-ET2890             |            |
| EC-720RI/L/M                                                                                                              | 12.80mm  | 3.8mm/1330 to 1990 | (8-14mm)            | 5mm 2400   | 727-ET2889             |            |
| EC-740T/M, L                                                                                                              | 10.7mm   | 3.2mm/1330 to 1690 | (8-14mm)            | 5mm 2400   | 727-ET2889             |            |
| EC-760P-V/M/L                                                                                                             | 11.1mm   | 3.2mm/1330 to 2010 | (8-14mm)            | 5mm 2400   | 727-ET2889             |            |
| EC-760R-V/I/L/M                                                                                                           | 12.00mm  | 3.8mm/1330 to 2010 | (8-14mm)            | 5mm 2400   | 727-ET2889             |            |
| EC-760S-A/M, L                                                                                                            | 12.8mm   | 3.8mm/1330 to 2010 | (8-14mm)            | 5mm 2400   | 727-ET2889             |            |
| EC-760S-G/L, M                                                                                                            | 12.8mm   | 3.8mm/1690-2010    | (8-14mm)            | 5mm 2400   | 727-ET2889             |            |
| EC-760Z-V/L/M                                                                                                             | 12.8mm   | 3.8mm/1330 to 2010 | (8-14mm)            | 5mm 2400   | 727-ET2889             |            |
| EC-760ZP-V/L/M                                                                                                            | 11.8mm   | 3.2mm/1330 to 2010 | (8-14mm)            | 5mm 2400   | 727-ET2889             |            |
| Dual Channel Scopes                                                                                                       |          |                    |                     |            |                        |            |
| EC-530DL/DM                                                                                                               | New      | 12.8mm             | 3.8/2.8mm 1690-2005 | (8-14mm)   | 5mm 2400               | TBC        |
| EG-530D                                                                                                                   | New      | 11.5mm             | 3.8/2.8mm 1100-1400 | (8-14mm)   | 5mm 2400               | TBC        |
| Duodenoscopes                                                                                                             |          |                    |                     |            |                        |            |
| ED-580T/580XT                                                                                                             |          | 13.1mm             | 4.2mm/1250-1550     | (8-14mm)   | 5mm 2400               | 727-ET2888 |
| Gastrosopes                                                                                                               |          |                    |                     |            |                        |            |
| EG-530CT                                                                                                                  |          | 10.8mm             | 3.8mm/1100-1400     | (8-14mm)   | 5mm 2400               | 727-ET2890 |
| EG-530FP                                                                                                                  |          | 8.5mm              | 2.8mm/1100-1400     | (8-14mm)   | 5mm 2400               | 727-ET2890 |
| EG-530N/NP                                                                                                                |          | 5.9mm              | 2.0mm/1100-1400     | (4-8mm)    | 3mm 1500               | 727-ET2490 |
| EG-530NW                                                                                                                  |          | 5.9mm              | 2.0mm/1100-1400     | (4-8mm)    | 3mm 1500               | 727-ET2490 |
| EG-530WR                                                                                                                  |          | 9.4mm              | 2.8mm/1100-1400     | (8-14mm)   | 5mm 2400               | 727-ET2890 |
| EG-580NW2                                                                                                                 |          | 5.9mm              | 2.4mm/1100-1400     | (4-8mm)    | 3mm 1500               | 727-ET2490 |
| EG-580RD                                                                                                                  |          | 9.8mm              | 3.2mm/1100-1400     | (8-14mm)   | 5mm 2400               | 727-ET2890 |
| EG-590WR                                                                                                                  |          | 9.6mm              | 2.8mm/1100-1400     | (8-14mm)   | 5mm 2400               | 727-ET2890 |
| EG-590ZW                                                                                                                  |          | 10.8mm             | 2.8mm/1100-1400     | (8-14mm)   | 5mm 2400               | 727-ET2890 |
| EG-600WI/L/M                                                                                                              |          | 12.0mm             | 3.8mm/1330-1990     | (8-14mm)   | 5mm 2400               | 727-ET2883 |
| EG-600WR                                                                                                                  |          | 9.3mm              | 2.8mm/1100-1400     | (8-14mm)   | 5mm 2400               | 727-ET2883 |
| EG-600ZW                                                                                                                  |          | 9.9mm              | 2.8mm/1100-1400     | (8-14mm)   | 5mm 2400               | 727-ET2883 |
| EG-720R                                                                                                                   |          | 9.3mm              | 2.8mm/1100-1400     | (8-14mm)   | 5mm 2400               | 727-ET2889 |
| EG-740N                                                                                                                   |          | 5.9mm              | 2.4mm/1100-1400     | (4-8mm)    | 3mm 1500               | 727-ET2489 |
| EG-760CT                                                                                                                  |          | 10.8mm             | 3.8mm/1100-1400     | (8-14mm)   | 5mm 2400               | 727-ET2889 |
| EG-760R                                                                                                                   |          | 9.3mm              | 2.8mm/1100-1400     | (8-14mm)   | 5mm 2400               | 727-ET2889 |
| EG-760Z                                                                                                                   |          | 9.9mm              | 2.8mm/1100-1400     | (8-14mm)   | 5mm 2400               | 727-ET2889 |

OTHER

Fujinon

**Sigmoidoscope**

|          |        |                |          |          |            |
|----------|--------|----------------|----------|----------|------------|
| ES-530WE | 12.8mm | 3.8mm/790-1090 | (8-14mm) | 5mm 2400 | 727-ET2890 |
|----------|--------|----------------|----------|----------|------------|

**Ultrasound Scopes**

|                 |        |                 |          |          |            |
|-----------------|--------|-----------------|----------|----------|------------|
| EB-530US Linear | 6.7mm  | 2.00mm/610-880  | (4-8mm)  | 3mm 1500 | 727-ET2493 |
| EG-580UT Linear | 13.9mm | 3.8mm/1250-1550 | (8-14mm) | 5mm 2400 | 727-ET2892 |
| EG-580UR        | 11.5mm | 2.8mm/1250-1550 | (8-14mm) | 5mm 2400 | 727-ET2891 |

Storz

**Bronchoscopes Video & Fibre**

Some scopes may require the cleaning adapter supplied with scopes to be used

|                 |       |               |         |            |             |
|-----------------|-------|---------------|---------|------------|-------------|
| 11301BNX        | 5.5mm | 2.3mm/650-930 | (4-8mm) | 3mm 1500   | 727-ET2485  |
| 11301 BNXK/BDXK | 5.5mm | 2.3mm/650-930 | (4-8mm) | 3mm 1500   | 727-ET2485  |
| 11302BNX/BDX    | 3.7mm | 1.5mm/650-940 | (4-8mm) | 1.5mm 1500 | 727-ET2485R |
| 11001BI/BN      | 5.2mm | 2.3mm/540-820 | (4-8mm) | 3mm 1500   | 727-ET2486  |

|                      |       |                  |         |            |             |
|----------------------|-------|------------------|---------|------------|-------------|
| 11002BD/11002BD2     | 3.6mm | 1.5mm/540-820    | (4-8mm) | 1.5mm 1500 | 727-ET2486R |
| 11003BC /11005BC     | 2.5mm | 1.2mm/540-820    | (1-4mm) | 1.5mm 1500 | 727-ET2186  |
| 11004BC/B1/11004BC1  | 6.4mm | 2.8mm/540-820    | (4-8mm) | 5mm 2400   | 727-ET2486L |
| 11009BC/BI           | 5.5mm | 2.6mm/540-820    | (4-8mm) | 5mm 2400   | 727-ET2486L |
| 11301AA/AA1/AB1/ABD1 | 2.8mm | 1.2mm/500 to 980 | (1-4mm) | 1.5mm 1500 | 727-ET2186  |
| 11301BN1 /BND1       | 5.2mm | 2.3mm/650-940    | (4-8mm) | 3mm 1500   | 727-ET2486  |
| 11302BD1/BD2/BDD2    | 3.7mm | 1.5mm/650-930    | (4-8mm) | 1.5mm 1500 | 727-ET2486R |

**Cystoscopes**

|               |        |           |         |            |            |
|---------------|--------|-----------|---------|------------|------------|
| 11272VP/VK/VN | 5.3mm  | 2.1mm/370 | (4-8mm) | 3mm 1500   | 727-ET2480 |
| 11272VH/VHK   | 5.2mm  | 2.1mm/370 | (4-8mm) | 3mm 1500   | 727-ET2485 |
| 11272C/CU     | 5.25mm | 2.1mm/370 | (4-8mm) | 3mm 1500   | 727-ET2480 |
| 11278AC/ACU   | 2.8mm  | 1.2mm/450 | (1-4mm) | 1.5mm 1500 | 727-ET2180 |

**Choledochoscopes**

|                     |       |           |         |            |            |
|---------------------|-------|-----------|---------|------------|------------|
| 11292 VSK/ADK1/DEK1 | 2.8mm | 1.2mm/500 | (1-4mm) | 1.5mm 1500 | 727-ET2180 |
| 11292 VP            | 5.3mm | 2.1mm/370 | (4-8mm) | 3mm 1500   | 727-ET2480 |

**Rectoscopes/Hysteroscopes**

|               |        |            |          |            |             |
|---------------|--------|------------|----------|------------|-------------|
| 13912 PKS/NKS | 11.8mm | 3.4mm/400  | (8-14mm) | 5mm 2400   | 727-ET2810  |
| 11264 BB1     | 3.5mm  | 1.48mm/240 | (4-8mm)  | 1.5mm 1500 | 727-ET2410R |

**Ureteroscopes**

|                |       |           |         |            |            |
|----------------|-------|-----------|---------|------------|------------|
| 11278A         | 2.2mm | 1.2mm/670 | (1-4mm) | 1.5mm 1500 | 727-ET2180 |
| 11278V/VK/ VSK | 2.8mm | 1.2mm/700 | (1-4mm) | 1.5mm 1500 | 727-ET2180 |

Wolf

**Bronchoscopes**

|         |       |             |         |            |            |
|---------|-------|-------------|---------|------------|------------|
| 7250001 | 5.2mm | 2.2mm/600mm | (4-8mm) | 3mm 1500   | 727-ET2482 |
| 7260001 | 6.2mm | 3.0mm/600mm | (4-8mm) | 3mm 1500   | 727-ET2482 |
| 7236001 | 2.9mm | 1.2mm/680mm | (1-4mm) | 1-5mm 1200 | 727-ET2182 |

**Cystoscope**

|          |       |             |         |          |            |
|----------|-------|-------------|---------|----------|------------|
| 7305.001 | 5.0mm | 2.5mm/400mm | (4-8mm) | 3mm 1500 | 727-ET2484 |
| 7305.006 | 5.0mm | 2.5mm/400mm | (4-8mm) | 3mm 1500 | 727-ET2484 |

**Cobra Vision**

|          |    |       |                   |         |            |            |
|----------|----|-------|-------------------|---------|------------|------------|
| 73561071 | !! | 3.3mm | 1.2mm&.08mm/680mm | (1-4mm) | 1.5mm 1500 | 727-ET2185 |
| 73561076 | !! | 3.3mm | 1.2mm&.08mm/680mm | (1-4mm) | 1.5mm 1500 | 727-ET2185 |

!! Do not supply a brush for 0.8mm, customer will need to use own