

Endotest Product Number Identification

The Endotest product codes have been set up dependant on the Endoscope Insertion tube tip size and Manufacturer type.

The Endotest Microbiological kit works the same way no matter what sample kit size you use, the only difference is your endoscope

To help understand the Microbiological sample kit you require, we have put together some simple explanations.

The following code depicts an Olympus Colonoscopes/Gastroscope Endotest Sample Kit code

727-ET2801

727

Just our bookkeeping code for finished goods products

ET

Endotest

28

V2 Version Pot Insertion tip size (used on scopes with IT sizes between 8-14mm)

- Col and gas scopes
- 24 = Scope insertion tube sizes between 4-8mm
- Bronchoscopes & Cystoscopes
- 21 = Scope insertion tube sizes between 1-4mm
- Small Diameter ureteroscopes etc

01

Scope Connection set pertaining to Olympus scopes
(01 to 49 = Olympus Kits)
Scope connection sets for Pentax scopes
(50 to 79 = Pentax Kits)
Scope connection sets for Storz, Fujinon, other
(80 to 90 = Pentax Kits)